



The device Audio Splint-Supratip is manufactured by Audio Technologies Srl  0373

The device Audio Kit-Supratip is assembled by Audio Technologies Srl



 audio<sup>®</sup>  
TECHNOLOGIES

**Audio Technologies S.r.l.**

Via Dell' Artigianato, 8 - 29020 Gossolengo, Piacenza, Italy

tel. +39.0523.778864 - fax +39.0523.770019

info@audiotechnologies.it - www.audiotechnologies.it



AUDIOSPLINT-*ST*

AUDIOKIT-*ST*



audio<sup>®</sup>  
TECHNOLOGIES

## Anti decubitus Self-adhesive Bonding System for Rhinoplasty

# AUDIOSPLINT-ST



### Audio Splint is composed by:

1. Anti-decubitus elastic polyurethane sheet. Self-adhesive surface to be applied over the first nasal dressing
2. Rigid element in light alloy. Self-adhesive surface to be applied over the soft element

### Soft element

Once the first nasal dressing is placed, put the self-adhesive side of the elastic sheet over the nasal dorsum and - while checking the profile - make sure the central flap of the splint adheres to the nasal dorsum.

Size	Packaging	Code
small	50 Soft elements + 50 Rigid elements	R 19.04 S ST
medium		R 19.05 M ST
large		R 19.06 L ST
small	10 Soft elements + 10 Rigid elements	R 19.04 S-P ST
medium		R 19.05 M-P ST
large		R 19.06 L-P ST

### Rigid element

The light-alloy splint keeps a pressure on the bony vault, easing the stabilization of a narrow nasal dorsum and a variable pressure on the cartilaginous vault.

The tip of the splint can adjust the pressure by keeping the nasal base in a narrow or wide position, at doctor's choice.

The central flap can be bent inwards to increase the supratip compression, avoiding the possibility to get edema and emphasizing the tip angle. If needed, the central strut can be trimmed using scissors.


### Use:

1. Place the first nasal medication
2. Place the self-adhesive soft element
3. Place the rigid element in aluminium

4. Place plasters and medications to keep the splint adherent to the nose.

**Removal:** Unstick the primary medication with a Joseph periosteal, paying attention not to elevate nasal skin.

**Post-surgery medication:** It's suggested to replace the splint every 5 - 7 days, for 3 weeks in total.



## Anti-decubitus Self-adhesive Bonding System for Rhynoplasty

# AUDIOKIT-ST

### Audio-Kit is composed by:

1. Soft polyurethane element with antidecubit cells. Self-adhesive surface to be applied over the sterile bandage
2. Rigid element in light alloy. Self-adhesive surface to be placed over the soft element
3. Sterile adhesive bandage for the first nasal dressing
4. Tincture of benzoin for increasing of the adhesion of the sterile adhesive bandage to the skin



Size	Packaging	Code
small	Each box contains 6 kits	R 19.01 K ST
medium	Each box contains 6 kits	R 19.02 K ST
large	Each box contains 6 kits	R 19.03 K ST

### Use:

1. Apply the first nasal dressing using steri-strips
2. Apply the anti-decubitus self-adhesive polyurethane soft element
3. Apply the self-adhesive aluminum rigid element
4. Apply overlapping strips to ensure the splint adheres on the nose

**Removal:** Use the blunt end of a Joseph-type elevator between the skin and the strips taking care not to detach the skin from the nasal framework.

**Post-operative follow up:** It is recommended to replace the dressing and the splint every 5 - 7 days, for a total of three medications.

