Audio Technologies srl
Via dell'Artigianato, 8
29020 Gossolengo
Piacenza - ITALY
www.audiotechnologies.it
info@audiotechnologies.it

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.....
The medical doctor is responsible for the correct procedures and surgical techniques. The instructions below are for information purposes only. The surgeon must evaluate the suitability of the procedures based on his training and experience.

AUDIO prostheses, for replacement of any or all of the bones of the ossicular chain, are manufactured from different materials: Platinum-Titanium-PTFE (fluoroplastic)-Porous polyethylene-Bioactive cera-

mic-Gold-Silicone-Nitinol. AUDIO prostheses have been successfully employed in case of otosclerosis, stapedo-vestibular ankylosis, outcomes of otitis media, trauma, middle ear congenital malformations. Possible complications due to the use of prostheses in the middle ear include: fixation, displacement, extrusion of the prosthesis, incus necrosis, perforation of the tympanic membrane, perilymph fistula, labyrinthitis, otitis media, occlusion of oval and round windows, formation of granuloma and allergic reaction to the material.

Attention: inform all patients who have had surgery which involved opening or removal of stapes footplate to avoid barotraumas which may cause penetration of the prosthesis into the oval window: inform patients who have had surgery not to touch the operated ear and not to introduce devices to clean the ear. Some models of prostheses are produced in various lengths, others must be trimmed during surgery. In all cases, the length of the prostheses must be carefully controlled with the measuring device (code 02.14). It is possible to reduce the risk of ejection of the device by inserting a sheet of autologous cartilage between the tympanum and the prosthesis.

Tympanic drainages - ATTENTION: if the tympanum is thin there is a risk of permanent perforation further to paracentesis and the risk that the device falls into the tympanic cavity: it is advisable to use suitable paracentesis instruments.

The tympanic drainage might be pushed in the tympanic cavity further to a trauma caused by ear cleaning devices.

Prostheses for stapedotomy - ATTENTION: Handle the prostheses with care, especially during trimming or when recovering the prosthesis from the tympanic cavity, avoiding any traction which could disconnect the PTFE stem from the prosthesis metal core.

Use the plate code SPL 03.06 to carry out the prosthesis trimming.

Prostheses with self-crimping hook: in some cases the ossicular chain is hypermobile. The pressure exerted by the operator on the incus when assembling the prosthesis might make the prosthesis stem sink into the internal ear. To avoid this problem assemble the hook of the self-crimping prosthesis to the incus before removing the stapes branches ("inversion of surgical steps") or hold up the long incus apophysis by means of a hook-shaped instrument (Audio Technologies ref. 02.11) when assembling the prosthesis hook. Use the Audio Technologies forceps ref. SPL 03.01 to assemble the prosthesis.

Cup-shaped prosthesis - ATTENTION - place a large strip of autologous tissue (vein or perichondrium) between the prosthesis and oval window. Ensure the prosthesis hook is securely anchored to the incus and there is no slipping of the prosthesis.

Elastic stem prostheses (TAP 07.18/PAP 07.19 - PAP 07.60 - TAP 07.71) - ATTENTION - having measured the distance to be protesised using measuring device 02.14, use an implant of equal length. Assemble the prosthesis very gently, avoiding traumatization of the stapes or relevant footplate. Implant the prosthesis fully extended. If the prosthesis is implanted in a compressed manner, it can exert excessive pressure on the stapes and cause dislocation or breakage.

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