


## Titanium and PTFE Ventilation Tubes

(PTFE: polytetrafluorethylene)

CODE TPL 04.02  
TPL 04.03

The edge of the external flange is easily gripped by microforceps during the insertion into the tympanic membrane. The large opening guarantees good drainage. Upon request, ventilation tubes are provided with a recovery thread.



TITANIUM	CODE (without recovery thread)	inner $\varnothing$ mm	flange $\varnothing$ mm	length mm	A mm
small	TPL 04.02.10 SF	1.0	2.25	1.55	1.15
medium	TPL 04.02.12 SF	1.25	2.50	1.55	1.15
large	TPL 04.02.15 SF	1.50	2.75	1.55	1.15

TITANIUM	CODE (with recovery thread)	inner $\varnothing$ mm	flange $\varnothing$ mm	length mm	A mm
small	TPL 04.02.10 CF	1.0	2.25	1.55	1.15
medium	TPL 04.02.12 CF	1.25	2.50	1.55	1.15
large	TPL 04.02.15 CF	1.50	2.75	1.55	1.15

PTFE-teflon	CODE (without recovery thread)	inner $\varnothing$ mm	flange $\varnothing$ mm	length mm	A mm
small	TPL 04.03.10 SF	1.0	2.25	1.55	1.15
medium	TPL 04.03.12 SF	1.25	2.50	1.55	1.15
large	TPL 04.03.15 SF	1.50	2.75	1.55	1.15

PTFE-teflon	CODE (with recovery thread)	inner $\varnothing$ mm	flange $\varnothing$ mm	length mm	A mm
small	TPL 04.03.10 CF	1.0	2.25	1.55	1.15
medium	TPL 04.03.12 CF	1.25	2.50	1.55	1.15
large	TPL 04.03.15 CF	1.50	2.75	1.55	1.15

Packaging: Dispenser box containing 6 sterile, individually packaged tubes.

